



A Predictive Diagnostic Test for Drug Treatment Outcomes in Advanced Breast Cancer Extended to Advanced Colon Cancer for Personalized Anticancer Therapy

Paul O.P. Ts'o PhD

Managing Director

CCC Diagnostics LLC
3918 Vero Rd. Suite B
Baltimore, MD 21227
www.cccdiag.com

PRESENTATIONS / ABSTRACTS



- K. H. Tkaczuk, N. S. Tait, O. Ioffe, M. Tan, M. Mohiuddin, S. Chumsri, D. A. VanEcho, M. J. Sutula, S. Lesko, S. Deamond, P. Ts'o; Drug Response Indicator Test (DRIT) as a predictive test for treatment outcomes in advanced breast cancer patients (ABC). Clin Oncol 27:15s, 2009 (suppl; abstr 1119) – ASCO 2009.
- N. B. Pandya, K. R. Tkaczuk, N. Tait, O. Ioffe, M. Tan, D. A. Van Echo, M. J. Sutula, S. A. Lesko, S. F. Deamond, P. O. Ts'o; A predictive test for therapeutic treatment outcomes of advanced gastrointestinal cancer patients (AGC). J Clin Oncol 27, 2009 (suppl; abstr e15079).

Personalized Anticancer Therapy

Background





- Cancer is a Highly Individualized Disease
- Breast Cancer = 184,450 new cases annually (2008)
40,930 deaths annually (2008)
- GI Cancer = 176,420 new cases annually
61,430 deaths annually
- Effective Treatment Relies on Drug Selection
 - Only 20-30% effective, significant side effects
- Some examples of Drug Response Indicator Guides include:
 - Estrogen receptor & antiestrogen therapy
 - Her2 & Trastuzumab therapy

Drug Response Indicator Test (DRIT)



- Diagnostic Service for Individual Patients
- Each DRI mechanistically related to specific drug
- Correlates DRI expression with clinical response
- Evaluates tumor response to 6 FDA approved anticancer drugs simultaneously
- Direct imaging of specified targets within a tumor tissue
- Selects effective drug treatments for Individual Patients
- Optimizes Use of Time Window for Treatment

DRIT - Methods



- Deparaffinized slides are subjected to antigen retrieval and stained with monoclonal antibodies.
- Antibodies directly labeled with Alexa dyes are used for staining
- Microscope calibrated to an external standard
- Stained slides are scanned to locate cytokeratin-positive cells.
- Five digital images are stored from different areas of each tumor. Also acquire images of the same fields of cells with filter cube C1 (DAPI) and C2 (cytokeratin).
- The fluorescence signals in each cropped image are analyzed by Image-Pro Plus software to select the intensities that separate the fluorescent objects (cells or cell clusters) from the background.
- IF expression in the tumor sections is measured in about 500 tumor cells utilizing 5 digital images of various areas of the tumor.
- Determinations of correlation are based on the relative level of comparative IF
- Results represented in Drug Response Indicator Units

Clinical Samples

- Clinical samples = slides cut from formaldehyde fixed, paraffin embedded breast cancer tissue
- Retrospective samples provided by:
 - University of Maryland Greenebaum Cancer Center , Baltimore, MD
 - Harborview Hospital , Baltimore , MD
- All clinical treatment outcomes determined via RECIST criteria



Drug Response Indicators Breast Cancer

- Drugs - FDA approved, NCCN recommended

DRUG	DRI
Trastuzumab	HER-2/neu
Antiestrogens	Estrogen Receptor
Taxanes/ Vinca Alkaloids	B-TUB III
5 FU	Thymidylate Synthase
Gemcitabine	Ribonucleotide Reductase

Drug Response Indicator Units



- DRI Molecular Unit - number of IgG-fluorescent conjugate molecules bound to the stained tissue that produces a fluorescence intensity of 2000
- Operational Unit - the value for one DRIU is 2000 F/P

Antibody Conjugate	# Conjugated antibody /pixel *
Anti- Beta-Tubulin (beta chain 4) terminus III – Alexa 532	3.3±0.09
Anti - Estrogen Receptor mouse clone TE111.5d11 – Alexa 532	10.2±0.25
Anti-thymidylate Synthase clone TS 106 – Alexa 532	8.8±0.27
Anti – Trastuzumab – Alexa 532	9.6±0.12

* number of antibody conjugates in a unit area of an image that yields a fluorescence intensity of one DRIU

DRI Expression Level Distribution

Breast Cancer



- Distribution of expression levels in sensitive vs resistant

TREATMENT	DRI	SENSITIVE	RESISTANT
		MEAN DRI UNITS	MEAN DRI UNITS
Hormonal	ER	> 9.7	< 9.7
Trastuzumab	HER-2/neu	> 7.5	< 7.5
Taxanes	TUB III	< 55	> 55
Capecitabine	TS	< 25	> 25

Eligibility Criteria – Breast Cancer



- Institutional IRB approved protocol/consent
- Stage 3 and 4 Breast Cancer Patients
- Mono or doublet therapy
- First two lines of therapy
- Neoadjuvant Therapy
- Excluded all experimental drug treatments
- Excluded all triplet therapy
- Excluded adjuvant therapy and 1 cycle treatments
- Provides basis for Prospective Study now underway at The University of Maryland Cancer Center

Patient Characteristics

Breast Cancer



	UMD	Harbor	Overall
<u>SEX RATIO</u> Male:Female	2/50	0/14	2/64
<u>RACE</u>			
African American	35	2	37
Caucasian	14	12	26
Asian	3	0	3
<u>AGE</u>			
Male	68 (58,78)		68 (58,78)
Female	54 (36,89)	61(41,89)	56 (36,89)
<u>STAGE – STUDY ENTRY</u>			
Stage III	1	4	5
Stage IV	51	10	61

DRIT Prediction Accuracy

Breast Cancer



UMD / Harbor Hospital Comparison

Treatment	U MD			Harbor Hospital		
	Sens. N=69	Resist	Predicted Response	Sens. N=22	Resist	Predicted Response
Hormonal Therapy	26/29	2/2	28/31	3/3	0/0	3/3
Mono- Chemotherapy	15/20	7/7	22/27	9/12	1/1	10/13
Doublet Chemotherapy	8/8	2/3	11/11	6/6	0/0	6/6
Overall Response	49/57	11/12	60/69	18/21	1/1	19/22
Prediction Accuracy %	86 %	92 %	87 %	86 %	100 %	86 %

Sens = sensitive ; Resis = Resistant

DRIT Prediction Accuracy

Breast Cancer

Hormonal Treatment



Treatment N=34	# Patients	DRIT PA	SOC PA	IHC PA
Aromatase Inhibitors	21	19/21	17/21	14/21
Tamoxifen	9	8/9	7/9	6/9
Faslodex	4	4/4	4/4	3/4
Total	34	31/34	28/34	23/34
% Response		91.2 %	82.4 %	67.6 %

PA= Predictive Accuracy; IHC = Immunohistochemistry; SOC = standard of care

DRIT Prediction Accuracy



Breast Cancer

Monochemotherapy/Trastuzumab

Treatment N=40	DRI	Prediction Accuracy (PA)		
		Sensitive Tumor	Resistant Tumor	Overall
Taxanes	Beta-Tubulin III	16/20	1/1	17/21
Capectabine	Thymidylate Synthase	7/11	3/3	10/14
Trastuzumab	HER-2/neu	1/1	4/4	5/5
Overall Response		24/32	8/8	32/40
% PA		75 %	100 %	80 %

Overall Results

Breast Cancer

DRIT vs Standard Care



Standard of Care (SOC) Treatment Success (TRS)					DRIT Prediction Accuracy (PA)			
Treatment	Resp	Non-Resp	Obj. Resp	% TRS	Sens	Resis	Predicted Response	PA %
Hormonal	28	6	28/34	82 %	29/32	2/2	31/34	90 %
Mono-Chemotherapy	24	16	24/40	56 %	24/32	8/8	32/40	82 %
Doublet Chemotherapy	14	3	14/17	73 %	14/14	2/3	16/17	91 %
Overall Response	66	25	66/91	73 %	67/78	12/13	79/91	87 %

DRIT Diagnostic Performance

Breast Cancer



N=91	Hormonal	Mono- Chemotherapy	Doublet Chemotherapy	Overall
Treatment interventions	34	40	17	91
Sensitivity	0.99	1	0.93	0.99
Specificity	0.4	0.50	1	0.52
Positive Predictive Value	0.86	0.74	1	0.86
Negative Predictive Value	0.92	1	0.67	0.92
Overall Accuracy	0.90	0.80	0.94	0.87

Conclusions



- DRIT measurement of estrogen receptor displays higher prediction accuracy than FDA approved commercial test
- DRIT displays a high degree of accuracy in predicting the response to non-hormonal therapies (chemotherapy & trastuzumab)
- DRIT predictions of treatment outcome outperformed the treatment success of the standard of care
- DRIT was able to identify ineffective treatments
- Preliminary results indicate that DRIT is applicable to a diverse population of patients



Drug Response Indicators GI Cancer

- Drugs - FDA approved, NCCN recommended

DRUG	DRI
Anthracycline	Topoisomerase II
Irinotecan	Topoisomerase I
Taxanes/ Vinca Alkaloids	B-TUB III
5 FU	Thymidylate Synthase
Gemcitabine	Ribonucleotide Reductase
Platinum Salts	ERCC1

DRIT Performance Accuracy

GI Cancer



Treatment Category	Sensitive	Resistant	Predicted Response	% Prediction Accuracy
Mono-chemotherapy	10/12	1/1	11/13	86.7 %
Doublet Chemotherapy	21/25	0/0	21/25	82.6 %
Triplet Chemotherapy	1/1	1/1	2/2	100 %
Overall Accuracy	32/38	2/2	34/40	85.0 %

DRIT Diagnostic Performance

GI Cancer



N=40	Mono- Chemotherapy	Doublet Chemotherapy	Triplet Chemotherapy	Overall
Treatment interventions	13	25	2	40
Sensitivity	1	1	1	1
Specificity	0.33	N/A	1	0.25
Positive Predictive Value	0.83	0.84	1	0.84
Negative Predictive Value	1	N/A	1	1
Overall Accuracy	0.87	0.83	1	0.85

Drug Response Indicator Test

Benefits for Patients



DRIT as a Diagnostic Service will:

- Individualize Anticancer Therapy Selection
- Lessen Exposure to Ineffective Anticancer Therapies
- Reduce the Need for Drug Combination Therapy
- Reduce Potential Side Effects to Patients
- Potential for Saving Significant Treatment Costs
- Improve Quality of Life

Drug Response Indicator Test

Cost Savings



- Based on Breast Cancer Retrospective Trials

Treatment	Resistant Tumors	Cost Savings
Endocrine (\$7800)	2	\$ 15,600
Taxane (\$12,000)	1	\$ 12,000
Capecitabine (\$ 15,000)	3	\$ 45,000
Herceptin (\$50,000)	4	\$ 200,000
Doublet Gemzar (\$20,000) Taxane / Herceptin (\$112,000)	2	\$ 132,000
Overall Savings (Drug Costs)	12	\$ 404,600
Ancillary Costs (\$4,000/3 cycles)		\$ 144,000
Total Savings		\$ 548,600
Savings/patient		\$ 8,000

Drug Response Indicator Test Cost Savings – Single Drug



- Saving Example based on single Drug = Taxotere
- Assumption = 50 % efficacy
- Drug Cost/ 100 patients = \$ 1,757,700
- Drug Savings / resistant patients = \$ 878,850
- Per patient Savings = \$ 8,788
- Savings/ US annual breast cancer population (40,000 patients)= \$ 351M

Drug Response Indicator Test

Competitive Advantage



- Rigorous in vitro qualification of all Drug Response Indicators
- Direct imaging of a specified target within a tissue sample
- Application to tissue sections prepared from paraffin blocks
- Standardization of fluorescent measurement
- Adaptable to new chemotherapeutic drugs
- Adaptable to additional cancers and biomarkers
- Predicts tumor response for up to six (6) FDA-Approved Drugs